

THAP12 RABBIT PAB

货号: S218810

产品全名: THAP12 兔多抗

基因符号 DAP4; THAP0; PRKRIR; P52rIPK

UNIPROT ID: O43422 (Gene Accession - BC021992)

背景: Upstream regulator of interferon-induced serine/threonine protein kinase R (PKR). May block the PKR-inhibitory function of DNAJC3, resulting in restoration of kinase activity and suppression of cell growth.

抗原: Fusion protein of human THAP12

经过测试的应用: ELISA, IHC

推荐稀释比: IHC: 40-200; ELISA: 5000-10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

亚型: Immunogen-specific rabbit IgG


纯化: Antigen affinity purification

种属反应性: Human, Mouse

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human ovarian cancer tissue using 218810(THAP12 Antibody) at a dilution of 1/45(Nucleus).