

## THIOREDOXIN 2 RABBIT PAB

货号: N225968

产品全名: Thioredoxin 2 兔多抗

基因符号 TXN; MTRX; TRX2; MT-TRX; COXPD29

**UNIPROT ID:** Q99757

**背景:** Thioredoxins (Trx) are small, multi-functional proteins with oxidoreductase activity and are ubiquitous in essentially all living cells. Trx contains a redox-active disulfide/dithiol group within the conserved Cys-Gly-Pro-Cys active site.

**抗原:** A synthesized peptide derived from human Thioredoxin 2

**经过测试的应用:** WB,IHC-P,ICC/IF,FC

**推荐稀释比:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 18 kDa; Observed MW: 12 kDa

**亚型:** IgG

**纯化:** Affinity Chromatography

**种属反应性:** Human, Mouse and Rat

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Signal Transduction

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Thioredoxin 2 in 3T3, K562 lysates using Thioredoxin 2 antibody.