

THIOREDOXIN DOMAIN CONTAINING PROTEIN 17 RABBIT PAB

货号: N225171

产品全名: Thioredoxin Domain Containing Protein 17 兔多抗

基因符号 TXNDC17; TXNL5; Thioredoxin domain-containing protein 17; 14 kDa thioredoxin-related protein; TRP14; Protein 42-9-9; Thioredoxin-like protein 5

UNIPROT ID: Q9BRA2

背景: Disulfide reductase. May participate in various redox reactions through the reversible oxidation of its active center dithiol to a disulfide and catalyze dithiol-disulfide exchange reactions. Modulates TNF-alpha signaling and NF-kappa-B activation. Has peroxidase activity and may contribute to the elimination of cellular hydrogen peroxide.

抗原: Synthesized peptide derived from TRP14 . at AA range: 20-100

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 14 kDa; Observed MW: 14 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Thioredoxin Domain Containing Protein 17 in HL-60, Raji, K562, SW480, A549, 22RV1 lysates using Thioredoxin Domain Containing Protein 17 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
