

THO COMPLEX SUBUNIT 1 RABBIT MAB

货号: N262627

产品全名: THO Complex Subunit 1 兔单克隆抗体

基因符号 Death domain containing protein p84N5; HPRI; hTRESK84;
Nuclear matrix protein p84; p84N5; THO complex 1; Tho1; Thoc1

UNIPROT ID: Q96FV9

背景: Regulates transcriptional elongation of a subset of genes.
Participates in an apoptotic pathway which is characterized by activation of caspase-6, increases in the expression of BAK1 and BCL2L1 and activation of NF-kappa-B.

抗原: A synthetic peptide of human Nuclear Matrix Protein p84

经过测试的应用: WB, ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R05-1A8

分子量: Calculated MW: 76 kDa; Observed MW: 84 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Tags & Cell Markers

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of Nuclear Matrix Protein p84 (green) in HeLa using Nuclear Matrix Protein p84 antibody, and DAPI (blue)



Western blot analysis of Nuclear Matrix Protein p84 in Jurkat, HeLa, COS7 lysates using THO Complex Subunit 1 antibody.