

## THROMBOSPONDIN 2 RABBIT PAB

货号: N225133

产品全名: Thrombospondin 2 兔多抗

基因符号 THBS2; TSP2; Thrombospondin-2

**UNIPROT ID:** P35442

背景: Adhesive glycoprotein that mediates cell-to-cell and cell-to-matrix interactions.

抗原: The antiserum was produced against synthesized peptide derived from the C-terminal region of human THBS2. AA range:1000-1050

经过测试的应用: WB,IHC-P,ELISA

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 130 kDa; Observed MW: 130 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded rat testis using Thrombospondin 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Thrombospondin 2 in mouse heart lysates using Thrombospondin 2 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---