

## THYMIDYLATE SYNTHASE RABBIT MAB

货号: N262994

产品全名: Thymidylate Synthase 兔单克隆抗体

基因符号 TS; TMS; HST422

**UNIPROT ID:** P04818

背景: Contributes to the de novo mitochondrial thymidylate biosynthesis pathway.

抗原: A synthetic peptide of human Thymidylate Synthase

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-2F2

分子量: Calculated MW: 36 kDa; Observed MW: 36 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human



成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling


储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Thymidylate Synthase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Thymidylate Synthase (green) in Jurkat using Thymidylate Synthase antibody, and DAPI (blue).



Western blot analysis of Thymidylate Synthase in K562 lysates using Thymidylate Synthase antibody.