

## TISSUE FACTOR RABBIT MAB

货号: N263581

产品全名: Tissue Factor 兔单克隆抗体

基因符号 F3; Tissue factor; TF; Coagulation factor III; Thromboplastin; CD antigen CD142

**UNIPROT ID:** P13726

**背景:** Initiates blood coagulation by forming a complex with circulating factor VII or VIIa. The [TF:VIIa] complex activates factors IX or X by specific limited proteolysis. TF plays a role in normal hemostasis by initiating the cell-surface assembly and propagation of the coagulation protease cascade.

**抗原:** Recombinant protein of human Tissue Factor

**经过测试的应用:** WB,IHC-P

**推荐稀释比:** WB: 1/500-1/1000 IHC: 1/50-1/100

**种属反应性:** Rabbit

**克隆性:** Rabbit Monoclonal

**克隆编号:** R07-3G5

**分子量:** Calculated MW: 33 kDa; Observed MW: 45 kDa

**亚型:** IgG


**纯化:** Affinity Purified

**种属反应性:** Human


**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Cardiovascular

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Tissue Factor antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Tissue Factor in A549, HL-60, U2OS lysates using Tissue Factor antibody.