

TISSUE FACTOR RABBIT MAB

货号: N263109

产品全名: Tissue Factor 兔单克隆抗体

基因符号 F3; Tissue factor; TF; Coagulation factor III; Thromboplastin; CD antigen CD142

UNIPROT ID: P13726

背景: Initiates blood coagulation by forming a complex with circulating factor VII or VIIa. The [TF:VIIa] complex activates factors IX or X by specific limited proteolysis. TF plays a role in normal hemostasis by initiating the cell-surface assembly and propagation of the coagulation protease cascade.

抗原: A synthetic peptide of human Tissue Factor

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-7G2

分子量: Calculated MW: 33 kDa; Observed MW: 45 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular


储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Tissue Factor antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Tissue Factor in mouse brain, rat brain lysates using Tissue Factor antibody.



Western blot analysis of Tissue Factor in HeLa, A549, MCF-7, U251, HL-60 lysates using Tissue Factor antibody.