

TNFAIP8 RABBIT MAB

货号: N261564

产品全名: TNFAIP8 兔单克隆抗体

基因符号 TNF alpha-induced protein 8; NDED; SCC-S2

UNIPROT ID: O95379

背景: Acts as a negative mediator of apoptosis and may play a role in tumor progression.

抗原: A synthetic peptide of human TNFAIP8

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-5F3

分子量: Calculated MW: 23 kDa; Observed MW: 21-23 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TNFAIP8 in CEM, K562 lysates using TNFAIP8 antibody.