

TNFAIP8L3 RABBIT PAB

货号: N225228

产品全名: TNFAIP8L3 兔多抗

基因符号 Tumor necrosis factor alpha-induced protein 8-like protein 3;
TNF alpha-induced protein 8-like protein 3; TNFAIP8-like protein 3

UNIPROT ID: Q5GJ75

背景: Acts as a lipid transfer protein.

抗原: Synthesized peptide derived from human TNFAIP8L3. at AA range:
15-45

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 33 kDa; Observed MW: 33 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat,Mouse

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol,
0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TNFAIP8L3 in 293T, mouse brain lysates using TNFAIP8L3 antibody.