

TOLL-LIKE RECEPTOR 4 RABBIT PAB

货号: N225529

产品全名: Toll-Like Receptor 4 兔多抗

基因符号 TLR4; Toll-like receptor 4; hToll; CD antigen CD284

UNIPROT ID: O00206

背景: Cooperates with LY96 and CD14 to mediate the innate immune response to bacterial lipopolysaccharide (LPS).

抗原: Recombinant Protein of TLR4

经过测试的应用: IHC-P

推荐稀释比: IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: -

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded mouse Spleen Tissue using Toll-Like Receptor 4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.