

TRANSCRIPTION INITIATION FACTOR IIB (2F6) MOUSE MAB

货号: N261397

产品全名: Transcription Initiation Factor IIB (2F6) 小鼠单抗

基因符号 General transcription factor IIB; General transcription factor TFIIB; gtf2b; RNA polymerase II transcription factor IIB; S300 II; S300-II; TF IIB; TF2B; TF2B_HUMAN; TFIIB; Transcription initiation factor IIB.

UNIPROT ID: Q00403

背景: General factor that plays a major role in the activation of eukaryotic genes transcribed by RNA polymerase II.

抗原: Purified recombinant human TFIIB protein fragments expressed in E.coli.

经过测试的应用: WB

推荐稀释比: WB: 1/500-1/1000

种属反应性: Mouse

克隆性: Mouse Monoclonal

克隆编号: 2F6-A3

分子量: Calculated MW: 35 kDa; Observed MW: 35 kDa

亚型: IgG2b


纯化: Affinity Purified

种属反应性: Human, Mouse, Monkey

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TFIIB in Jurkat, MCF-7, COS7, 3T3 and HeLa lysates using TFIIB antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
