

## TRANSCRIPTION INITIATION FACTOR IIB RABBIT MAB

货号: N263576

产品全名: Transcription Initiation Factor IIB 兔单克隆抗体

基因符号 General transcription factor IIB; General transcription factor TFIIB; gtf2b; RNA polymerase II transcription factor IIB; S300 II; S300-II; TF IIB; TF2B; TF2B\_HUMAN; TFIIB; Transcription initiation factor IIB.

**UNIPROT ID:** Q00403

背景: General factor that plays a major role in the activation of eukaryotic genes transcribed by RNA polymerase II.

抗原: A synthetic peptide of human TFIIB

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R06-4C2

分子量: Calculated MW: 35 kDa; Observed MW: 35 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics&Histone Acetylation

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TFIIIB in HeLa, A549, HL-60, U2OS, C6 lysates using TFIIIB antibody.



Western blot analysis of TFIIIB in 3T3, Raw264.7 lysates using Transcription Initiation Factor IIB antibody.