

TRANSFERRIN RABBIT MAB

货号: N263015

产品全名: Transferrin 兔单克隆抗体

基因符号 TF; Serotransferrin; Transferrin; Beta-1 metal-binding globulin; Siderophilin

UNIPROT ID: P02787

背景: Transferrins are iron binding transport proteins which can bind two Fe³⁺ ions in association with the binding of an anion, usually bicarbonate. It is responsible for the transport of iron from sites of absorption and heme degradation to those of storage and utilization. Serum transferrin may also have a further role in stimulating cell proliferation.

抗原: A synthetic peptide of human Transferrin

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R03-4K6

分子量: Calculated MW: 77 kDa; Observed MW: 77 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Transferrin in K562 lysates using Transferrin antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
