

TRANSFERRIN RECEPTOR 1 RABBIT MAB

货号: N263016

产品全名: Transferrin Receptor 1 兔单克隆抗体

基因符号 TFRC; Transferrin receptor protein 1; TR; TfR; TfR1; Trfr; T9; p90; CD71

UNIPROT ID: P02786

背景: Binds the telomeric double-stranded TTAGGG repeat and negatively regulates telomere length. Involved in the regulation of the mitotic spindle. Component of the shelterin complex (telosome) that is involved in the regulation of telomere length and protection.

抗原: Recombinant protein of human Transferrin Receptor

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-6H2

分子量: Calculated MW: 85 kDa; Observed MW: 90 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of Transferrin Receptor (green) in HeLa using Transferrin Receptor antibody, and DAPI (blue)



Western blot analysis of Transferrin Receptor in K562, HeLa lysates using Transferrin Receptor 1 antibody.