

TRANSGLUTAMINASE 2 RABBIT MAB

货号: N263017

产品全名: Transglutaminase 2 兔单克隆抗体

基因符号 TGM2; Protein-glutamine gamma-glutamyltransferase 2;
Tissue transglutaminase; Transglutaminase C; TG(C); TGC; TGase C;
Transglutaminase H; TGase H; Transglutaminase-2; TGase-2

UNIPROT ID: P21980

背景: Catalyzes the cross-linking of proteins and the conjugation of polyamines to proteins.

抗原: A synthetic peptide of human Transglutaminase 2

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-8K1

分子量: Calculated MW: 77 kDa; Observed MW: 77 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology


储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Transglutaminase 2 in rat Brain lysates using Transglutaminase 2 antibody



Immunohistochemistry analysis of paraffin-embedded mouse colon using Transglutaminase 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Transglutaminase 2 in HUVEC lysates using Transglutaminase 2 antibody.