

TRANSGLUTAMINASE 2 RABBIT PAB

货号: N225824

产品全名: Transglutaminase 2 兔多抗

基因符号 TGM2; Protein-glutamine gamma-glutamyltransferase 2; Tissue transglutaminase; Transglutaminase C; TG(C); TGC; TGase C; Transglutaminase H; TGase H; Transglutaminase-2; TGase-2

UNIPROT ID: P21980

背景: Catalyzes the cross-linking of proteins and the conjugation of polyamines to proteins.

抗原: A synthesized peptide derived from human TGM2

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 77 kDa; Observed MW: 77 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human,Mouse

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TGM2 in HUVEC lysates using Transglutaminase 2 antibody.