

TRANSMEMBRANE PROTEIN 72 RABBIT PAB

货号: N225244

产品全名: Transmembrane Protein 72 兔多抗

基因符号 Transmembrane protein 72; Kidney-specific secretory protein of 37 kDa; C10orf127; KSP37

UNIPROT ID: A0PK05

背景: This gene encodes a transmembrane protein which may be expressed specifically in the kidney.

抗原: Synthesized peptide derived from human Tmem72 Polyclonal

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 30 kDa; Observed MW: 30 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Transmembrane Protein 72 in various lysates using Transmembrane Protein 72 antibody.