

TRANSMEMBRANE PROTEIN 85 RABBIT MAB

货号: N263004

产品全名: Transmembrane Protein 85 兔单克隆抗体

基因符号 PIG17; TMEM85

UNIPROT ID: Q5J8M3

背景: May mediate anti-apoptotic activity.

抗原: A synthetic peptide of human TMEM85

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-6D5

分子量: Calculated MW: 20 kDa; Observed MW: 20 kDa

亚型: IgG

纯化: Affinity Purified

种属反应性: Human, Mouse and Rat,Hamster


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TMEM85 in Jurkat, rat Brain, C6, CHO-K1, HeLa



3T3 lysates using Transmembrane Protein 85 antibody.