

TRIMETHYL-HISTONE H4 (LYS20) RABBIT PAB

货号: N225636

产品全名: TriMethyl-Histone H4 (Lys20) 兔多抗

基因符号 H4K20me3; H4F2; HIST; HIST1H4A; Histone H4

UNIPROT ID: P62805

背景: The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a member of the histone H4 family.

抗原: A synthesized peptide derived from human Histone H4 (tri methyl K20)

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 11 kDa; Observed MW: 11 kDa

亚型: IgG

纯化: Affinity Chromatography


种属反应性: Human, Mouse and Rat

Modification: Methylated


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human stomach using Histone H4 (tri methyl K20) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Histone H4 (tri methyl K20) in HeLa lysates using TriMethyl-Histone H4 (Lys20) antibody.