

TRMT11 RABBIT PAB

货号: N225202

产品全名: TRMT11 兔多抗

基因符号 TRMT11; C6orf75; MDS024; tRNA; guanine(10)-N2)-methyltransferase homolog; tRNA guanosine-2'-O-methyltransferase TRMT11 homolog

UNIPROT ID: Q7Z4G4

背景: Catalytic subunit of an S-adenosyl-L-methionine-dependent tRNA methyltransferase complex that mediates the methylation of the guanosine nucleotide at position 10 (m²G₁₀) in tRNAs.

抗原: Synthesized peptide derived from the C-terminal region of human TRMT11.

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 53 kDa; Observed MW: 53 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: -

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TRMT11 in HepG2, JK, Colo, HeLa lysates using TRMT11 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
