

TRPV5 RABBIT MAB

货号: N263030

产品全名: TrpV5 兔单克隆抗体

基因符号 CAT2; ECAC1; OTRPC3

UNIPROT ID: Q9NQA5

背景: This gene is a member of the transient receptor family and the TrpV subfamily. The calcium-selective channel encoded by this gene has 6 transmembrane-spanning domains, multiple potential phosphorylation sites, an N-linked glycosylation site, and 5 ANK repeats. This protein forms homotetramers or heterotetramers and is activated by a low internal calcium level.

抗原: A synthetic peptide of human TRPV5

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R03-7A4

分子量: Calculated MW: 83 kDa; Observed MW: 83 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience&Transient Receptor Potential Channel

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TRPV5 in K562, C6, 3T3, HeLa lysates using TrpV5 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
