

## TRYPTASE DELTA RABBIT PAB

货号: N225019

产品全名: Tryptase delta 兔多抗

基因符号 Tryptase delta (EC 3.4.21.59; Delta-tryptase; HmMCP-3-like tryptase III; Mast cell mMCP-7-like; Tryptase-3)

**UNIPROT ID:** Q9BZJ3

背景: Tryptase is the major neutral protease present in mast cells and is secreted upon the coupled activation-degranulation response of this cell type.

抗原: Synthesized peptide derived from Tryptase-3 at AA range: 51-100

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 27 kDa; Observed MW: 27 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: -

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Tryptase delta in HEPG2 lysates using Tryptase delta antibody.