

TSG101 RABBIT MAB

货号: N263032

产品全名: TSG101 兔单克隆抗体

基因符号 TSG101; Tumor susceptibility gene 101 protein; ESCRT-I complex subunit TSG101

UNIPROT ID: Q99816

背景: Component of the ESCRT-I complex, a regulator of vesicular trafficking process. Binds to ubiquitinated cargo proteins and is required for the sorting of endocytic ubiquitinated cargos into multivesicular bodies (MVBs). Mediates the association between the ESCRT-0 and ESCRT-I complex. Required for completion of cytokinesis; the function requires CEP55.

抗原: A synthetic peptide of human TSG101

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-9A7

分子量: Calculated MW: 44 kDa; Observed MW: 44 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TSG101 in K562 lysates using TSG101 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
