

## TSG101 RABBIT PAB

货号: N225751

产品全名: TSG101 兔多抗

基因符号 TSG101; Tumor susceptibility gene 101 protein; ESCRT-I complex subunit TSG101

**UNIPROT ID:** Q99816

**背景:** Component of the ESCRT-I complex, a regulator of vesicular trafficking process. Binds to ubiquitinated cargo proteins and is required for the sorting of endocytic ubiquitinated cargos into multivesicular bodies (MVBs). Mediates the association between the ESCRT-0 and ESCRT-I complex. Required for completion of cytokinesis; the function requires CEP55.

**抗原:** A synthesized peptide derived from human TSG101

**经过测试的应用:** WB,IHC-P,ICC/IF,FC

**推荐稀释比:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 44 kDa; Observed MW: 44 kDa

**亚型:** IgG


**纯化:** Affinity Chromatography

**种属反应性:** Human, Mouse and Rat

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Cell Biology

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TSG101 in Jurkat lysates using TSG101 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---