

## UCHL5 RABBIT MAB

货号: N263049

产品全名: UCHL5 兔单克隆抗体

基因符号 UCH37; CGI-70; INO80R; UCH-L5

**UNIPROT ID:** Q9Y5K5

背景: Protease that specifically cleaves 'Lys-48'-linked polyubiquitin chains.

抗原: A synthetic peptide of human UCH37

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R03-9C1

分子量: Calculated MW: 38 kDa; Observed MW: 38 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology


储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of UCH37 (green) in Hela using UCH37 rat spleen lysates using UCHL5 antibody, and DAPI (blue)



Western blot analysis of UCH37 in rat spleen lysates using UCHL5 antibody.



Western blot analysis of UCH37 in 3T3, Hela lysates using UCH37 antibody.