

UFC1 RABBIT MAB

货号: N263052

产品全名: UFC1 兔单克隆抗体

基因符号 NEDSG; HSPC155

UNIPROT ID: Q9Y3C8

背景: E2-like enzyme which forms an intermediate with UFM1 via a thioester linkage.

抗原: Recombinant protein of human UFC1

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R06-5E3

分子量: Calculated MW: 19 kDa; Observed MW: 19 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of UFC1 (green) in HeLa using UFC1 antibody, and DAPI (blue)



Western blot analysis of UFC1 in K562, HeLa lysates using UFC1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
