

UFC1 RABBIT MAB

货号: N263051

产品全名: UFC1 兔单克隆抗体

基因符号 NEDSG; HSPC155

UNIPROT ID: Q9Y3C8

背景: E2-like enzyme which forms an intermediate with UFM1 via a thioester linkage.

抗原: Recombinant protein of human UFC1

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-5E1

分子量: Calculated MW: 19 kDa; Observed MW: 19 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using UFC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of UFC1 in K562, rat Brain, C6, 3T3 lysates using UFC1 antibody.