

## USO1 RABBIT MAB

货号: N263068

产品全名: USO1 兔单克隆抗体

基因符号 TAP; VDP; P115

**UNIPROT ID:** O60763

**背景:** General vesicular transport factor required for intercisternal transport in the Golgi stack; it is required for transcytotic fusion and/or subsequent binding of the vesicles to the target membrane. May well act as a vesicular anchor by interacting with the target membrane and holding the vesicular and target membranes in proximity.

**抗原:** A synthetic peptide of human Vesicle docking protein p115

**经过测试的应用:** WB,IP

**推荐稀释比:** WB: 1/500-1/1000 IP: 1/20

**种属反应性:** Rabbit

**克隆性:** Rabbit Monoclonal

**克隆编号:** R07-4B8

**分子量:** Calculated MW: 108 kDa; Observed MW: 108 kDa

**亚型:** IgG


**纯化:** Affinity Purified

**种属反应性:** Mouse

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Signal Transduction

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of USO1 in 3T3 lysates using USO1 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---