

USP13 RABBIT MAB

货号: N263057

产品全名: USP13 兔单克隆抗体

基因符号 Deubiquitinating enzyme 13; Isopeptidase T3; ISOT3; UBPI3; USP13

UNIPROT ID: Q92995

背景: Mediates deubiquitination of BECN1, a key regulator of autophagy, leading to stabilize the PIK3C3/VPS34-containing complexes. Also deubiquitinates USP10, an essential regulator of p53/TP53 stability.

抗原: A synthetic peptide of human USP13

经过测试的应用: WB,ICC/IF,IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-8D3

分子量: Calculated MW: 97 kDa; Observed MW: 97 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of USP13 in U2OS lysates using USP13 antibody.



Immunocytochemistry analysis of USP13 (green) in HeLa using USP13 antibody, and DAPI (blue).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
