

VASP RABBIT MAB

货号: N263061

产品全名: VASP 兔单克隆抗体

基因符号 Vasodilator-stimulated phospho蛋白

UNIPROT ID: P50552

背景: VASP (vasodilator-stimulated phosphoprotein), is involved in the maintenance of cytoarchitecture by interacting with actin-like filaments. VASP shares a limited degree of homology with the amino terminus of WASP, which is frequently mutated in WAS patients. An established substrate of cAMP and cGMP dependent kinases, VASP is phosphorylated on a regulatory serine residue 157 and localizes to focal adhesions, microfilaments and highly active regions of the plasma membrane.

抗原: A synthetic peptide of human VASP

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-2B0

分子量: Calculated MW: 40 kDa; Observed MW: 46 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of VASP in HeLa lysates using VASP antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
