

## VAV2 RABBIT MAB

货号: N263604

产品全名: VAV2 兔单克隆抗体

基因符号 VAV2; Vav 2 oncogene; VAV-2

**UNIPROT ID:** P52735

背景: Vav proteins belong to the Dbl family of guanine nucleotide exchange factors (GEFs) for Rho/Rac small GTPases. The three identified mammalian Vav proteins (Vav1, Vav2 and Vav3) differ in their expression. Vav1 is expressed only in hematopoietic cells and is involved in the formation of the immune synapse. Vav2 and Vav3 are more ubiquitously expressed.

抗原: A synthetic peptide of human VAV2

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R06-8F5

分子量: Calculated MW: 101 kDa; Observed MW: 101 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of VAV2 in HeLa, MCF-7, HepG2, U2OS, U251 lysates using VAV2 antibody.



Immunocytochemistry analysis of VAV2 (green) in HeLa using VAV2 antibody, and DAPI (blue).