

VCAM1 RABBIT MAB

货号: N261757

产品全名: VCAM1 兔单克隆抗体

基因符号 VCAM1; LICAM; VCAM 1; V CAM 1; DKFZp779G2333; CD106

UNIPROT ID: P29533

背景: Important in cell-cell recognition. Appears to function in leukocyte-endothelial cell adhesion. Interacts with integrin alpha-4/beta-1 (ITGA4/ITGB1) on leukocytes, and mediates both adhesion and signal transduction. The VCAM1/ITGA4/ITGB1 interaction may play a pathophysiologic role both in immune responses and in leukocyte emigration to sites of inflammation.

抗原: Recombinant protein of mouse VCAM1

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-6E9

分子量: Calculated MW: 81 kDa; Observed MW: 100 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Mouse

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of VCAM1 in 3T3 lysates using VCAM1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
