

VITAMIN D BINDING PROTEIN RABBIT MAB

货号: N263074

产品全名: Vitamin D Binding Protein 兔单克隆抗体

基因符号 GC; Vitamin D-binding protein; DBP; VDB; Gc-globulin; Group-specific component

UNIPROT ID: P02774

背景: DBP associates with membrane-bound immunoglobulin on the surface of B-lymphocytes and with IgG Fc receptor on the membranes of T-lymphocytes.

抗原: A synthetic peptide of human Vitamin D Binding 蛋白

经过测试的应用: WB, ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R03-5G8

分子量: Calculated MW: 53 kDa; Observed MW: 53 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Vitamin D Binding protein in rat Brain, 3T3, HeLa lysates using Vitamin D Binding protein antibody.



Immunocytochemistry analysis of Vitamin D Binding Protein (green) in HEPG2 using Vitamin D Binding Protein antibody, and DAPI (blue).