

## WDR5 RABBIT PAB

货号: N225917

产品全名: WDR5 兔多抗

基因符号 WDR5; BIG3; WD repeat-containing protein 5; BMP2-induced 3-kb gene 蛋白

**UNIPROT ID:** P61964

**背景:** Contributes to histone modification. May position the N-terminus of histone H3 for efficient trimethylation at 'Lys-4'. As part of the MLL1/MLL complex it is involved in methylation and dimethylation at 'Lys-4' of histone H3. H3 'Lys-4' methylation represents a specific tag for epigenetic transcriptional activation.

**抗原:** A synthesized peptide derived from human WDR5

**经过测试的应用:** WB,IHC-P,FC

**推荐稀释比:** WB: 1/500-1/1000 IHC: 1/50-1/100 FC: 1/50-1/100

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 37 kDa; Observed MW: 37 kDa

**亚型:** IgG


**纯化:** Affinity Chromatography

**种属反应性:** Human, Mouse and Rat

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Epigenetics and Nuclear Signaling

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of WDR5 in HeLa lysates using WDR5 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---