

## WNT16 RABBIT MAB

货号: N263083

产品全名: Wnt16 兔单克隆抗体

基因符号 WNT16; Protein Wnt-16

**UNIPROT ID:** Q9UBV4

背景: Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters .

抗原: A synthetic peptide of human Wnt16

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-2G2

分子量: Calculated MW: 41 kDa; Observed MW: 41 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Stem Cells

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Wnt16 in K562, rat Brain lysates using Wnt16 antibody.



Western blot analysis of Wnt16 in mouse brain lysates using Wnt16 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---