

WNT2B RABBIT PAB

货号: N225941

产品全名: Wnt2b 兔多抗

基因符号 WNT13; Wnt2b; XWNT2

UNIPROT ID: Q93097

背景: Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters.

抗原: A synthesized peptide derived from human Wnt2b

经过测试的应用: WB,ICC/IF,FC

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 44 kDa; Observed MW: 44 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Stem Cells

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Wnt2b in HepG2 lysates using Wnt2b antibody.