

XCT RABBIT MAB

货号: N263087

产品全名: xCT 兔单克隆抗体

基因符号 CCBRI; Cysteine/glutamate transporter; SLC7A11; solute carrier family 7; xCT

UNIPROT ID: Q9UPY5

背景: Sodium-independent, high-affinity exchange of anionic amino acids with high specificity for anionic form of cystine and glutamate.

抗原: A synthetic peptide of human xCT

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-2C5

分子量: Calculated MW: 55 kDa; Observed MW: 55 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Mouse

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of xCT in K562 lysates using xCT antibody.