

YY1 RABBIT MAB

货号: N263098

产品全名: YY1 兔单克隆抗体

基因符号 YY1; Delta transcription factor; INO80 complex subunit S; NF-E1; Yin and yang 1; YY-1

UNIPROT ID: P25490

背景: Its activity is regulated by transcription factors and cytoplasmic proteins that have been shown to abrogate or completely inhibit YY1-mediated activation or repression. For example, it acts as a repressor in absence of adenovirus E1A protein but as an activator in its presence. Acts synergistically with the SMAD1 and SMAD4 in bone morphogenetic protein (BMP)-mediated cardiac-specific gene expression.

抗原: A synthetic peptide of human YY1

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-1D3

分子量: Calculated MW: 45 kDa; Observed MW: 65 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of YY1 in C6, 3T3, HeLa lysates using YY1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
