

## YY1 RABBIT MAB

货号: N263097

产品全名: YY1 兔单克隆抗体

基因符号 YY1; Delta transcription factor; INO80 complex subunit S; NF-E1; Yin and yang 1; YY-1

**UNIPROT ID:** P25490

**背景:** Its activity is regulated by transcription factors and cytoplasmic proteins that have been shown to abrogate or completely inhibit YY1-mediated activation or repression. For example, it acts as a repressor in absence of adenovirus E1A protein but as an activator in its presence. Acts synergistically with the SMAD1 and SMAD4 in bone morphogenetic protein (BMP)-mediated cardiac-specific gene expression.

**抗原:** A synthetic peptide of human YY1

**经过测试的应用:** WB,IP

**推荐稀释比:** WB: 1/500-1/1000 IP: 1/20

**种属反应性:** Rabbit

**克隆性:** Rabbit Monoclonal

**克隆编号:** R02-8A3

**分子量:** Calculated MW: 45 kDa; Observed MW: 65 kDa

**亚型:** IgG


**纯化:** Affinity Purified

**种属反应性:** Human, Mouse and Rat

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Epigenetics and Nuclear Signaling

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of YY1 in K562, rat Brain, C6, 3T3, HeLa lysates using YY1 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---