

HIS MAB

产品名称: His

货号: 26904

规格: 100 μ L

Symbol: His-tag

描述: HIS Tag 小鼠单抗

背景: The His-tag antibody recognizes any protein containing continuous 6 or more histidine residues.

Immunogen: His synthetic peptide

经过测试的应用: ELISA, WB, IF, IHC

推荐稀释比: ELISA 1:2000-1:5000 WB 1:1000-1:5000 IF 1:100-1:200

Concentration: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG


纯化: Purified from ascites

Preservative: No

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Recognizes any protein containing continuous 6 histidine residues.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing



蛋白质印迹:

蛋白质印迹: 重组的HIS标签融合蛋白, 用HIS鼠单克隆抗体(1:1000)稀释, 蛋白上样量为50ng进行蛋白印迹分析