

ARF6 MAB

产品名称: Arf6 小鼠单抗

货号: 26029

基因符号: ARF6

描述: Arf6 小鼠单抗

背景: ARF6 (ADP-ribosylation factor 6) is a member of the ARF super-family. ARF genes encode small GTPases that increase the ADP-ribosyltransferase activity of cholera toxin and are critical for vesicular trafficking as activators of phospholipase D. Arf6 regulates membrane trafficking and the actin cytoskeleton at the plasma membrane and functions as a regulatory molecule of phagocytosis.

免疫原: Recombinant full length protein of Arf6.

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释比: ELISA: 1:100-1:5000WB: 1:100-1:1000IHC: 1:10-1:1000

Concentration: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal


亚型: IgG

纯化: Purified from ascites

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Arf6 antibody recognizes Arf6 of vertebrates.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing



免疫沉淀/蛋白质印迹:

was used for western blot (Cat# 26029).

HEK-293 whole-cell lysate

