

CD26 MAB (BEGELOMAB)

货号: 28146

产品全名: CD26 (Begelomab) 单抗

基因符号 ADABP;ADCP2;CD26;DPPIV;TPI03

描述: CD26 (begelomab) 单抗

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Mus musculus

亚型: IgG2b

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 CD26

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8 % trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. CD26 protein is highly expressed on the surface of Expi293 cell membrane. Flow cytometry analysis with 1 $\mu\text{g}/\text{mL}$ Anti-CD26 (begelomab biosimilar) mAb (28146) on Expi293 cells transfected with Human CD26 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram), and Isotype antibody on Expi293 transfected with irrelevant protein (Orange histogram)



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
