

BCMA MAB (BIOSIMILAR MAB)

货号: 28025

产品全名: BCMA (Biosimilar) 单抗

基因符号 CD269;TNFRSF17;BCM;BCMA

描述: BCMA () 单抗

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: ELISA; Flow Cyt

推荐稀释比: ELISA 1:5000-10000; Flow Cyt 1:100

种属反应性: Humanized

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 BCMA

成分: Lyophilized from sterile PBS, pH 7.4. 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. ELISA plate pre-coated by 2µg/ml (100µl/well) Human BCMA, hFc-His Tag 11129 can bind Anti-BCMA (Neutralizing antibody clone huC11D5.3) in a linear range of 3.71-22.29 ng/ml.




Figure 2. Flow cytometry analysis with 15 µg/mL Anti-BCMA mAb (28025) on Expi293 cells transfected with Human BCMA protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).