

BIOTINYLATED ANTI-GPRC5D (DM61) RABBIT MAB

货号: 28366

产品全名: Biotinylated Anti-GPRC5D(DM61) 兔单克隆抗体

基因符号 GPRC5D

描述: Biotinylated Anti-GPRC5D antibody(DM61) 兔单克隆抗体

背景: The protein encoded by this gene is a member of the G protein-coupled receptor family; however, the specific function of this gene has not yet been determined.

经过测试的应用: ELISA; Flow Cyt

推荐稀释比: IHC 1:200, ELISA 1:5000-10000; Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 GPRC5D

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).