

BIOTINYLATED ANTI-PSCA (DM87) RABBIT MAB

货号: 28416

产品全名: Biotinylated Anti-PSCA(DM87) 兔单克隆抗体

基因符号 PSCA;UNQ206;PRO232

描述: Biotinylated Anti-PSCA antibody(DM87) 兔单克隆抗体

背景: This gene encodes a glycosylphosphatidylinositol-anchored cell membrane glycoprotein. In addition to being highly expressed in the prostate it is also expressed in the bladder; placenta; colon; kidney; and stomach. This gene is up-regulated in a large proportion of prostate cancers and is also detected in cancers of the bladder and pancreas. This gene includes a polymorphism that results in an upstream start codon in some individuals; this polymorphism is thought to be associated with a risk for certain gastric and bladder cancers. Alternative splicing results in multiple transcript variants.

经过测试的应用: ELISA; Flow Cyt

推荐稀释比: ELISA 1:5000-10000; Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 PSCA

成分: Lyophilized from sterile PBS, pH 7.4. 5 % - 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. 1×10^5 of Expi 293 cell line transfected with human PSCA were stained with 100 μ l of 1:100 diluted Biotinylated Anti-PSCA antibody (DM87), Rabbit mAb (Red histogram) or isotype control (Blue histogram) respectively, washed and then stained with Streptavidin APC. The experimental samples were analyzed by flow cytometry.