

PCSK9 MAB (BOCOCIZUMAB)

货号: 28096

产品全名: PCSK9(Bococizumab) 单抗

基因符号 FH3;FHCL3;HCHOLA3;LDLCQ1;NARC-1;NARC1;PC9

描述: PCSK9(bococizumab) 单抗

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Humanized

亚型: IgG2

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 PCSK9

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis under cell membrane permeable condition with 1 $\mu\text{g}/\text{mL}$ Anti-PCSK9 (bococizumab biosimilar) mAb (28096) on Expi293 cells transfected with Human PCSK9 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).