

CANIS LUPUS FAMILIARIS ANXA2 PROTEIN, HIS TAG

货号: 12124

产品全名: Canis Lupus Familiaris ANXA2 蛋白

规格: 10/50/100 µg

基因符号 Annexin A2;Annexin-2

目标蛋白: ANXA2

UNIPROT ID: Q6TEQ7

描述: Recombinant canis ANXA2 protein with C-terminal 6xHis tag

背景: This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. This protein functions as an autocrine factor which heightens osteoclast formation and bone resorption. This gene has three pseudogenes located on chromosomes 4, 9 and 10, respectively. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. Annexin A2 expression has been found to correlate with resistance to treatment against various cancer forms. [provided by RefSeq, Dec 2019]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 39.3 kDa after removal of the signal peptide. The apparent molecular mass of dANXA2-His is approximately 35-55 kDa due to glycosylation.

分子特征: ANXA2(Ser2-Asp339) 6xHis tag

纯化: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Canis ANXA2 Protein, His Tag on SDS-PAGE under reducing condition.