

CD96 (DMC285) IGG1 CHIMERIC MAB

货号: 28193

产品全名: CD96(DMC285) IgG1 Chimeric 单抗

基因符号 TACTILE

描述: CD96 antibody(DMC285) IgG1 Chimeric 单抗

背景: The protein encoded by this gene belongs to the immunoglobulin superfamily. It is a type I membrane protein. The protein may play a role in the adhesive interactions of activated T and NK cells during the late phase of the immune response. It may also function in antigen presentation. Alternative splicing generates multiple transcript variants encoding distinct isoforms.

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit:Human Fc chimeric IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 CD96

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with Anti-CD96 (DMC285) on Expi293 cells transfected with human CD96 (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).