

## CD70 MAB (CHONGQING PRECISION MAB)

货号: 28159

产品全名: CD70 (Chongqing Precision) 单抗

基因符号 CD27L;CD27LG;TNFSF7

描述: CD70 (Chongqing Precision) 单抗

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Humanized

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 CD70

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at  $-80^{\circ}\text{C}$  (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1 $\mu\text{g}/\text{mL}$  Anti-CD70 (Chongqing Precision biosimilar) mAb (28159) on Expi293 cells transfected with Human CD70 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---